

---

**Product Data Sheet**

---

Product Name: Signal Enhancer Hikari

Cat. No.: GC14265

**Chemical Properties**

Cas. No.

Formula  $C_{13}H_{14}N_4OS$ 

M.Wt 274.34

Solubility Soluble to 20 mM in DMSO and to 10 mM in ethanol Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure 

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA