
Product Data Sheet

Product Name: RD3-0028

Cat. No.: GC32340

Chemical Properties

Cas. No. 3886-39-3

SMILES C12=CC=CC=C1CSSC2Formula C8H8S2

M.Wt

168.28

Solubility Soluble in DMSO

Storage

Store at -20°C

General tips

For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Protocol**

Cell experiment: 50 mL of HeLa cell suspensions (30000 cells/well) and the Long strain of RSV(25 TCID50/well) in MEM supplemented with 0.1% bovine serum albumin and antibiotics are added to each well in a 96-well round-bottomed microtiter plate that is filled with 50mL of MM in the presence or absence of several concentrations of RD3-0028. The anti-RSV assay is performed primarily with the MTT method[1].

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

Product Data Sheet

Animal experiment:

Mouse: Aerosols are generated from reservoirs containing 0.3 to 7.0 mg of RD3-0028 per mL. Solutions of ribavirin are prepared in saline containing 2.5 to 60 mg/mL. Mice are treated intraperitoneally with 100 mg of cyclophosphamide per kg of body weight 5 days before virus inoculation. The mice are weighed, anesthetized with sodium pentobarbital (50 mg/kg), and inoculated intranasally with approximately 10⁵ PFU of RSV A2 in 50 µL (day 0). From day 1 through day 3, the mice are exposed to the RD3-0028 or ribavirin aerosol. Placebo consisted of 10% DMSO-saline containing 1% Tween 80. On day 4, the day on which untreated mice had the maximum RSV pulmonary titer, all animals are killed and the lungs of each mouse are removed and for virus quantification[2].

Rat: 14C-RD3-0028 is dissolved in 10% dimethyl sulphoxide (DMSO)/saline containing 1% Tween 80 for administration by the aerosol route. The aerosol is generated using a head-exposure chamber, mono-position, with a mist generator. The mass median aerodynamic diameter of the aerosol particle is 2.1 µm. Rats are exposed to the 14C-RD3-0028 aerosol for 15min and killed at indicated times after the end of exposure. Because the aerosol treatment with 14C-RD3-0028 exposed the rat to 8.8 mg (91.5 kBq)/animal, the solution of this compound is orally administered at a dose of 8.8 mg/body[3].

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

Product Data Sheet

References:

- [1]. Watanabe W, et al. Novel anti-respiratory syncytial(RS) viral compounds: benzodithiin derivatives. Biochem Biophys Res Commun. 1998 Aug 28;249(3):922-6.
- [2]. Sudo K, et al. Efficacy of RD3-0028 aerosol treatment against respiratory syncytial virus infection in immunosuppressed mice. Antimicrob Agents Chemother. 1999 Apr;43(4):752-7.
- [3]. Sudo K, et al. Pharmacokinetics of a benzodithiin (RD3-0028) following aerosol treatment in rat.

Background

RD3-0028 is a potent and selective inhibitor of RSV replication with an EC50 of 4.5 μ M.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

Product Data Sheet

RD3-0028 has a 50% effective concentration of 4.5 μM and a 50% cytotoxic concentration of 271.0 μM which is superior to that of ribavirin. RD3-0028 inhibits different RSV strains at a low concentration (4.5-11.0 μM) using the MTT method. Using the MTT method, EC50 values of RD3-0028 against tested strains are lower than those of ribavirin. RD3-0028 does not inhibit the replication of measles virus, influenza A virus, herpes simplex virus types 1 and 2, or human cytomegalovirus[1].

Aerosols generated from reservoirs containing RD3-0028 (7 mg/mL) administered for 2 h twice daily for 3 days significantly reduces the pulmonary titer of RSV-infected mice. It is clear that the minimal effective dose of RD3-0028 for RSV-infected mice is significantly less than that of ribavirin, the only compound currently available for use against RSV disease. Furthermore, the RD3-0028 aerosol administration protect the lungs of infected, CYP-treated mice against tissue damage, as evidenced by the preservation of the lung architecture and a reduction in pulmonary inflammatory infiltrates. RD3-0028 aerosol is not toxic for mice at the therapeutic dose[2]. The plasma concentration of RD3-0028 is maintained at the same level from 5 min to 1 h, and decreases with a half-life of 2.2 h for 1 ± 8 h. The excretion of radioactivity in the urine and faeces at 24 h after aerosol treatment is 89.3 and 4.5%, respectively, indicating that almost all the radioactivity is rapidly excreted in the urine. The excretion of total radioactivity is 98.9% within 168 h[3].

[1]. Watanabe W, et al. Novel anti-respiratory syncytial(RS) viral compounds: benzodithiin derivatives. *Biochem Biophys Res Commun*. 1998 Aug 28;249(3):922-6. [2]. Sudo K, et al. Efficacy of RD3-0028 aerosol treatment against respiratory syncytial virus infection in immunosuppressed mice. *Antimicrob Agents Chemother*. 1999 Apr;43(4):752-7. [3]. Sudo K, et al. Pharmacokinetics of a benzodithiin (RD3-0028) following aerosol treatment in rat.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA
