
Product Data Sheet

Product Name: RNase Inhibitor, Murine

Cat. No.: GK12001

Components

Components	Size
RNase Inhibitor, Murine(40 U/ul)	3000 U (75 ul)
RNase Inhibitor, Murine(40 U/ul)	15000 U (375 ul)

Store the components at -20°C.

Features

Applications	The product inhibits RNases A, B and C specifically. Improved resistance to oxidation compared to human/porcine RNase inhibitor. Ideal for reactions when low DTT concentration is necessary (e.g., Real-time PCR).
Shipping	Ship with blue ice.
Storage	Store the components at -20°C.
Usage	For Research Use Only! Not For Use in Humans.

Protocol

1. Used in IVT(In vitro translation), RT-PCR, cDNA synthesis and other application Add RNase Inhibitor, Murine to achieve a final concentration of 1U/μl in the reaction. Be careful that RNase Inhibitor should be added before other components that may possibly containing RNase (i.e. enzymes, plasmid) Murine RNase Inhibitor is a 50 kDa protein that will be inactivated in denaturing conditions such as high temperature and denaturing agents, it must be used below 50°C.

2. Used in adding poly(A) tail to RNA In the experiment of adding poly(A) tail to RNA, RNase Inhibitor can be added to enhance the stability of the RNA. 0.5 μl of RNase inhibitor to a 20 μl reaction is appropriate when the reaction is set-up. The additional volume can be subtracted from the amount of H2O used in the reaction.

Background

RNase Inhibitor, Murine is a 50 kD recombinant protein isolated from E. coli containing the ribonuclease inhibitor gene originally from mouse. Our product RNase inhibitor has specificity for RNases A, B and C, but couldn't inhibit RNase 1, RNase T1, S1 Nuclease,

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

Product Data Sheet

RNase H or RNase from Aspergillus. It can bind several kinds of RNase non-covalently in a 1:1 ratio with high affinity. Especially, according to previous researches, no inhibition for enzyme activity is observed when RNase Inhibitor is used within AMV or M-MuLV reverse transcriptases, taq DNA polymerase or SP6/ T7/ T3 RNA polymerases, means it can be taken in various experiments to prevent RNA degradation.

Compared to RNase from human, RNase inhibitor from mouse lacks the pair of cysteines which is very sensitive to oxidation and lead to inactivation of the inhibitor. As a consequence, RNase Inhibitor, Murine significantly improves its resistance to oxidation compared to RNase inhibitors from human/porcine. It is stable at low DTT concentrations (less than 1 mM), and this character makes it preferred for reactions where low concentration DTT is necessary to the reaction (eg. Real-time RT-PCR).

Our product can be used in many kinds of reactions to prevent RNA degradation, for example, RT-PCR, cDNA synthesis, In vitro transcription, enzymatic RNA labeling reaction or any other applications.

RNase Inhibitor, Murine 50 kD
RNases A, B, C, RNase 1, RNase T1, S1, RNase H, Aspergillus
RNase, RNA, 1:1, AMV, M-MuLV, taq
DNA, SP6/T7/T3 RNA, RNA
RNase Inhibitor
RNase, DTT, 1 mM, DTT, RT-PCR
RNA, RT-PCR, cDNA, RNA

Caution: Product has not been fully validated for medical applications. For research use only.
Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com
Address: 10292 Central Ave. #205, Montclair, CA, USA