

---

## Product Data Sheet

---

Product Name: Dengue Envelope-2 32kDa

Cat. No.: GP25530

Batch No.: 1

### Product Data

Purity	>98%	Source	E.coli.
Physical Appearance	solid	Shipping Condition	Shipped with Ice Packs.
Formulation	Phosphate buffer.		

### Introduction

Dengue fever is affected by 1 of 4 closely linked virus serotypes of the genus Flavivirus, family Flaviviridae. One might have dengue fever infected by the different serotype virus after the primary infection. Detection of particular antibodies to dengue viruses is used for the diagnosis in clinic. Recently, lateral flow rapid test products have become a most suitable and known method in the clinical diagnosis. Though, there is difficulty for manufacturers to have the dengue antigens with a whole coverage for Dengue IgG & IgM recognition for all 4 serotype infections as well as with colloid gold binding ability. 8 dengue antigens have been established for the lateral flow products which are well defined for their coverage of dengue IgG & IgM recognition. The researcher can select the products below based on their own specific application.

### Stability

Dengue Envelope-2 32kDa although stable at 4°C for 1 week, should be stored below -18°C .Please prevent freeze thaw cycles.

### Background

Recombinant Dengue Virus Subtype 2 Envelope 32kDa, produced in E.coli, amino acids 45-297 and fused to a 6xHis Tag, containing domains I+II of the dengue envelope.

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA