

---

**Product Data Sheet**

---

Product Name: H-Phe(4-I)-OH

Cat. No.: GA10046

**Chemical Properties**

Cas. No. 24250-85-9

Chemical Name (2S)-2-amino-3-(4-iodophenyl)propanoic acid

SMILES C1=CC(=CC=C1CC(C(=O)O)N)IFormula  $C_9H_{10}INO_2$ 

M.Wt 291.1

Solubility  $H_2O$  : 20 mg/mL (68.71 mM; ultrasonic and adjust pH to 2 with HCl)

Storage

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

H-Phe(4-I)-OH

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA