
Product Data Sheet

Product Name: HMP Linker

Cat. No.: GA10056

Chemical Properties

Cas. No. 68858-21-9

Chemical Name 2-[4-(hydroxymethyl)phenoxy]acetic acid

SMILES C1=CC(=CC=C1CO)OCC(=O)OFormula $C_9H_{10}O_4$

M.Wt 182.2

Solubility $\geq 18.2\text{mg/mL}$ in DMSO, $\geq 30.75\text{ mg/mL}$ in EtOH

Storage

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

HMP Linker

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA