

---

**Product Data Sheet**

---

Product Name: Fmoc-Asp(OtBu)-Wang resin

Cat. No.: GA10064

**Chemical Properties**

Cas. No. N/A

Chemical Name (S,E)-tert-butyl 3-((((9H-fluoren-9-yl)methoxy)(hydroxy)methylene)amino)-4-oxobutanoate

SMILES CC(C)(OC(C[C@@])(/N=C(OCC1C2=CC=CC=C2C3=CC=CC=C31)\O)([H])C=O)=O)CFormula C<sub>23</sub>H<sub>25</sub>NO<sub>5</sub> M.Wt 395.45

Solubility Storage

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

Fmoc-Asp(OtBu)-Wang resin

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA