
Product Data Sheet

Product Name: Fmoc-D-Ile-Wang resin

Cat. No.: GA10258

Chemical Properties

Cas. No. N/A

Chemical Name (R,E)-(9H-fluoren-9-yl)methyl hydrogen (1-oxohexan-2-yl)carbonimidate

SMILES CCCC[C@]1(/N=C(OCC1C2=CC=CC=C2C3=CC=CC=C31)\O)([H])C=OFormula $C_{21}H_{23}NO_3$ M.Wt 337.41

Solubility Storage

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

Fmoc-D-Ile-Wang resin

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA