
Product Data Sheet

Product Name: Fmoc-Arg(Pbf)-Wang resin

Cat. No.: GA10269

Chemical Properties

Cas. No. N/A

Chemical Name 9H-fluoren-9-ylmethyl N-[(2S)-5-[[amino-[(2,2,4,6,7-pentamethyl-3H-1-benzofuran-5-yl)sulfonylamino]methylidene]amino]-1-oxopentan-2-yl]carbamate

SMILES CC1=C2C(=C(C(=C1C)S(=O)C(=O)NC(=NCCCC(C=O)NC(=O)OCC3C4=CC=CC=C4C5=CC=CC=C35)N)C)CC(O2)(C)C

Formula M.Wt

Solubility Storage Store at -20°C

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT, or blue ice upon request.

Structure

Background

Fmoc-Arg(Pbf)-Wang resin

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA
