
Product Data Sheet

Product Name: HMPA-AM Resin

Cat. No.: GA10421

Chemical Properties

Cas. No. N/A

Formula 100-200:mesh,1%DVB

M.Wt 0.5~1.5mmol/g

Solubility

Storage

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

HMPA-AM Resin

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA