
Product Data Sheet

Product Name: H-Trp-2-Chlorotrityl Resin

Cat. No.: GA10688

Chemical Properties

Cas. No. N/A

Chemical Name (2S)-(2-chlorophenyl)(phenyl)(p-tolyl)methyl 2-amino-3-(1H-indol-3-yl)propanoate

SMILES CC(C=C1)=CC=C1C(C2=CC=CC=C2)(OC([C@H](CC3=CNC4=CC=CC=C34)N)=O)C5=C(Cl)C=CC=C5

Formula 100-200:mesh,1%DVB M.Wt 0.3-0.8mmol/g

Solubility Storage

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

Background

H-Trp-2-Chlorotrityl Resin

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA