

---

## Product Data Sheet

---

Product Name: Fmoc-His(Trt)-Wang resin

Cat. No.: GA10702

### Chemical Properties

Cas. No. N/A

Chemical Name (S)-(9H-fluoren-9-yl)methyl (1-oxo-3-(1-trityl-1H-imidazol-4-yl)propan-2-yl)carbamate

SMILES  
O=C(N([H])[C@@](C([H])=O)([H])C([H])  
([H])C(N=C1[H])=C([H])N1C(C2=C([H])C([H])=C([H])C([H])=C2[H])  
(C3=C([H])C([H])=C([H])C([H])=C3[H])C4=C([H])C([H])=C([H])C([H])=C4[H])OC([H])  
([H])C5([H])C6=C([H])C([H])=C([H])C([H])=C6C7=C([H])C([H])=C([H])C([H])=C57

Formula M.Wt

Solubility Storage 2-8°C, Sealed storage.

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

### Background

Fmoc-His(Trt)-Wang resin

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA