
Product Data Sheet

Product Name: Fmoc-Trp-Wang resin

Cat. No.: GA10851

Chemical Properties

Cas. No. N/A

Chemical Name 9H-fluoren-9-ylmethyl N-[(2S)-1-(1H-indol-3-yl)-3-oxopropan-2-yl]carbamate

SMILES C1=CC=C2C(=C1)C(C3=CC=CC=C32)COC(=O)NC(CC4=CNC5=CC=CC=C54)C=O

Formula 100-200:mesh,1%DVB M.Wt 0.3-0.8mmol/g

Solubility Storage

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

Background

Fmoc-Trp-Wang resin

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA