
Product Data Sheet

Product Name: H-Phe-2-Chlorotrityl Resin

Cat. No.: GA11166

Chemical Properties

Cas. No. N/A

Chemical Name (S)-2-amino-3-phenylpropanoic acid

SMILES O=C([C@@](C([H])([H])C1=C([H])C([H])=C([H])C([H])=C1[H])([H])N([H])[H])O[H]

Formula 100-200:mesh,1%DVB M.Wt 0.3-0.8mmol/g

Solubility Storage

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

H-Phe-2-Chlorotrityl Resin

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA