
Product Data Sheet

Product Name: Fmoc-Argininol(Tos)

Cat. No.: GA11169

Chemical Properties

Cas. No.

Chemical Name (S)-(9H-fluoren-9-yl)methyl (1-hydroxy-5-(3-tosylguanidino)pentan-2-yl)carbamate

SMILES O=C(N[C@H](CO)CCCNC(NS(=O)(=O)C1=CC=C(C)C=C1)=O)=N)OCC2C3=CC=CC=C3C4=CC=CC=C42Formula $C_{28}H_{32}N_4O_5S$ M.Wt 536.0

Solubility Soluble in water or 1% acetic acid Storage Desiccate at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

Fmoc-Argininol(Tos)

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA