

---

**Product Data Sheet**

---

Product Name: Fmoc-Dap(Dnp)-OH

Cat. No.: GA11215

**Chemical Properties**

Cas. No. 140430-54-2

Chemical Name (2S)-3-(2,4-dinitroanilino)-2-(9H-fluoren-9-ylmethoxycarbonylamino)propanoic acid

SMILES C1=CC=C2C(=C1)C(C3=CC=CC=C32)COC(=O)NC(CNC4=C(C=C(C=C4)[N+](=O)[O-])[N+](=O)[O-])C(=O)OFormula C<sub>24</sub>H<sub>20</sub>N<sub>4</sub>O<sub>8</sub> M.Wt 492.44

Solubility Storage

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

Fmoc-Dap(Dnp)-OH

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA