

---

**Product Data Sheet**

---

Product Name: Wang Resin

Cat. No.: GA11357

**Chemical Properties**

Cas. No. 201058-08-4

Chemical Name [4-(Hydroxymethyl)phenoxyethyl]polystyrene

SMILES OCC1=CC=C(OCC2=CC=C(C=C2)[\*])C=C1

Formula 100-200:mesh,1%DVB M.Wt 1.0~1.2mmol/g

Solubility Storage

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

Wang Resin

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA