
Product Data Sheet

Product Name: ([ring-D₅]Phe⁶)-Somatostatin-14

Cat. No.: GA20016

Chemical Properties

Cas. No. 1926163-72-5

Formula C₇₆H₉₉D₅N₁₈O₁₉S₂

M.Wt 1642.94

Solubility Soluble in DMSO

Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

Stable isotope-labeled Somatostatin-14 useful for pharmacokinetic/pharmacodynamic studies in combination with glpbio product Somatostatin-14 or the corresponding generic API (Somatostatin Ph. Eur.).

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA