
Product Data Sheet

Product Name: (Arg⁸)-Vasopressin (4-9)

Cat. No.: GA20033

Chemical Properties

Cas. No. 96027-30-4

Formula C₂₅H₄₄N₁₂O₈S

M.Wt

672.77

Solubility Soluble in DMSO

Storage

Store at -20°C

General
tips

For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

Main metabolite of (Arg⁸)-Vasopressin. It improved the scopolamine-induced deficit of spatial memory in rats about 10'000 fold more than (Arg⁸)-Vasopressin and therefore has the most potent activity of all the endogenous metabolites of (Arg⁸)-Vasopressin

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA