
Product Data Sheet

Product Name: (D-Arg⁰,Hyp³,β-(2-thienyl)-Ala^{5·8},D-Phe⁷)-Bradykinin

Cat. No.: GA20071

Chemical Properties

Cas. No. 103412-36-8

Formula C₅₆H₈₃N₁₉O₁₃S₂

M.Wt

1294.53

Solubility Soluble in DMSO

Storage

Store at -20°C

General tips

For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

D-Arg-Arg-Pro-Hyp-Gly-Thi-Ser-D-Phe-Thi-Arg, a potent antagonist of bradykinin. Inhibitor of human urinary kallikrein, Ki 0.5 μM.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA