
Product Data Sheet

Product Name: (D-Arg²)-Kyotorphin

Cat. No.: GA20076

Chemical Properties

Cas. No. 70904-57-3

Formula C₁₅H₂₃N₅O₄ M.Wt 337.38

Solubility Soluble in DMSO Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

D-Kyotorphin showed higher analgesic activity in rats than kyotorphin, though the analgesic actions of these diastereomeric dipeptides probably resulted from different mechanisms. Contrary to kyotorphin, only minimal modulation of the Na/Cl dependant opi

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA