
Product Data Sheet

Product Name: (D-Cys(tBu)²,Thr(tBu)⁶)-Leu-Enkephalin-Thr

Cat. No.: GA20081

Chemical Properties

Cas. No. 131766-24-0

Formula C₄₁H₆₂N₆O₉S M.Wt 815.04

Solubility Soluble in DMSO Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

The tBu-modified hexapeptide BUBUC is the most potent δ -receptor-specific enkephalin analog.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA