

---

**Product Data Sheet**

---

Product Name: (Deamino-Cys<sup>3</sup>,Nle<sup>4</sup>,Arg<sup>5</sup>,D-2-Nal<sup>7</sup>,Cys<sup>11</sup>)- $\alpha$ -MSH (3-11) amide

Cat. No.: GA20085

**Chemical Properties**

Cas. No. 613222-50-7

Formula C<sub>56</sub>H<sub>76</sub>N<sub>18</sub>O<sub>9</sub>S<sub>2</sub> M.Wt 1209.47

Solubility Soluble in DMSO Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

JKC-366, potent analog of HS024, is a selective MC4 receptor antagonist.

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA