
Product Data Sheet

Product Name: (D-His²)-LHRH

Cat. No.: GA20134

Chemical Properties

Cas. No. 53634-19-8

Formula C₅₅H₇₅N₁₇O₁₃ M.Wt 1182.31

Solubility Soluble in DMSO Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure (D-His²)-LHRH**Background**

GnRH agonist with residual LH-releasing activity (10 %) on rat pituitary cells. Possible contaminant of synthetic gonadorelin due to the high propensity of His derivatives to racemize during coupling.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA