
Product Data Sheet

Product Name: (H-Asp-Glu-Val-Asp)₂-Rhodamine 110

Cat. No.: GA20202

Chemical Properties

Cas. No. 220846-75-3

Formula C₅₆H₆₆N₁₀O₂₃ M.Wt 1247.19

Solubility Soluble in DMSO Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

(DEVD)₂-R110, non-fluorescent peptidyl (DEVD)-Rhodamine 110-bisamide, is a highly sensitive, photostable fluorogenic substrate for the determination of caspase-3 activity. Characteristic for substrates based on Rhodamine 110 (R110) is that they contain

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA