
Product Data Sheet

Product Name: (Ile⁷⁶)-TNF- α (70-80) (human)

Cat. No.: GA20226

Chemical Properties

Cas. No. 163045-82-7

Formula C₅₅H₉₁N₁₅O₁₆ M.Wt 1218.42

Solubility Soluble in DMSO Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure (Ile⁷⁶)-TNF- α (70-80) (human)**Background**

The peptide PSTHVLITHTI corresponds to residues 70-80 of human TNF- α substituted with Ile at position 76 and exhibits many of the biological activities found for TNF- α . (Ile⁷⁶)-TNF- α (70-80) has been shown to enhance human polymorphonucle

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA