

---

**Product Data Sheet**

---

Product Name: (Leu<sup>31</sup>,Pro<sup>34</sup>)-Peptide YY (human)

Cat. No.: GA20230

**Chemical Properties**

Cas. No. 179986-95-9

Formula C<sub>195</sub>H<sub>296</sub>N<sub>54</sub>O<sub>56</sub>

M.Wt

4292.83

Solubility Soluble in DMSO

Storage

Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

LP-NYY, a Y1-subtype-selective agonist of the NPY receptor.

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA