
Product Data Sheet

Product Name: (Leu⁸,D-Trp²²,Tyr²⁵)-Somatostatin-28

Cat. No.: GA20231

Chemical Properties

Cas. No. 77909-99-0

Formula C₁₃₈H₂₀₉N₄₁O₄₀S₂

M.Wt 3146.56

Solubility Soluble in DMSO

Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

This somatostatin analog which can be iodinated has been used to characterize somatostatin binding sites in rat brain and hamster insulinoma membranes.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA