
Product Data Sheet

Product Name: (Lys(Ac)^{5.8.12.16})-Histone H4 (1-25)-Gly-Ser-Gly-Ser-Lys(biotinyl)

Cat. No.: GA20233

Chemical Properties

Cas. No. 2022956-65-4

Formula C₁₄₁H₂₄₃N₅₃O₄₃S

M.Wt

3400.88

Solubility Soluble in DMSO

Storage

Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA