
Product Data Sheet

Product Name: (Lys(Me)₃²⁷)-Histone H3 (21-44)-Gly-Lys(biotinyl)

Cat. No.: GA20238

Chemical Properties

Cas. No. 2022956-67-6

Formula C₁₃₀H₂₂₁N₄₄O₃₃S

M.Wt

2960.53

Solubility Soluble in DMSO

Storage

Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA