
Product Data Sheet

Product Name: (Met¹⁸⁶)-Melanocyte Protein PMEL 17 (185-193) (human, bovine, mouse)

Cat. No.: GA20255

Chemical Properties

Cas. No. 181477-43-0

Formula C₄₇H₇₄N₁₀O₁₄S M.Wt 1035.23

Solubility Soluble in DMSO Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

Background

This modified sequence of the gp100 epitope gp(209-217), in which a methionine replaced the natural threonine at position 2, bound to the HLA-A2 molecule with greater affinity than the unmodified peptide and was shown to have an increased ability to ge

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA