
Product Data Sheet

Product Name: (Phenylac¹,D-Tyr(Et)²,Lys⁶,Arg⁸,des-Gly⁹)-Vasopressin

Cat. No.: GA20274

Chemical Properties

Cas. No. 129520-65-6

Formula C₅₄H₇₆N₁₄O₁₁ M.Wt 1097.29

Solubility Soluble in DMSO Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

Selective and highly potent V receptor antagonist (antivasopressor) with an antiantidiuretic/antivasopressor ED ratio of about 570.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA