
Product Data Sheet

Product Name: (Propionyl¹,D-Tyr(Et)²,Val⁴,Abu⁶,Arg^{8,9})-Vasopressin

Cat. No.: GA20281

Chemical Properties

Cas. No. 121250-95-1

Formula C₅₃H₈₂N₁₆O₁₁ M.Wt 1119.34

Solubility Soluble in DMSO Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

This linear vasopressin analog is as potent a V (antidiuretic) antagonist as its cyclic parent compound. Its V (vasopressor) antagonism, however, is significantly reduced, resulting in a promising anti-V/anti-V selectivity.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA