
Product Data Sheet

Product Name: (Tyr⁰)-C-Type Natriuretic Peptide (32-53) (human, porcine, rat)

Cat. No.: GA20323

Chemical Properties

Cas. No. 142878-79-3

Formula C₁₀₂H₁₆₆N₂₈O₃₀S₃

M.Wt 2360.81

Solubility Soluble in DMSO

Storage Store at -20°C

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT, or blue ice upon request.

Structure

Background

This analog of C-Type Natriuretic Peptide (32-53) (human, porcine, rat) was synthesized for iodination.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA