

---

**Product Data Sheet**

---

Product Name: (Tyr<sup>1</sup>)-TRAP-7

Cat. No.: GA20327

**Chemical Properties**

Cas. No. 149440-16-4

Formula C<sub>45</sub>H<sub>67</sub>N<sub>11</sub>O<sub>10</sub> M.Wt 922.1

Solubility Soluble in DMSO Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

The TRAP-7 analog YFLLRNP is an antagonist to  $\alpha$ -thrombin and to the thrombin receptor agonist peptide SFLLRNP in human platelets. It might be a useful tool for the differentiation between several possible activation states of the human platelet th

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA