
Product Data Sheet

Product Name: (Tyr³⁶)-pTH-Related Protein (1-36) (human, mouse, rat)

Cat. No.: GA20331

Chemical Properties

Cas. No. 213779-11-4

Formula C₁₉₄H₃₀₃N₅₉O₅₃

M.Wt 4309.91

Solubility DMSO : 100 mg/mL (23.20 mM; Need ultrasonic) Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

Background

(Tyr³)-pTHrP (1-36) has been used for radioiodination.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA