
Product Data Sheet

Product Name: 3-(4-Methoxy-benzylsulfanyl)-3-methyl-butyric acid

Cat. No.: GA20392

Chemical Properties

Cas. No. 268219-99-4

Formula $C_{13}H_{18}O_3S$ M.Wt 254.35

Solubility Soluble in DMSO Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

S-protected derivative for the introduction of a deamino-penicillamine residue. Building block for the preparation of oxytocin and vasopressin antagonists. Peptide thioesters for chemical ligation may be obtained by Fmoc-SPPS using a deaminopenicillami

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA