
Product Data Sheet

Product Name: 4-(Fmoc-hydrazino)-benzoic acid

Cat. No.: GA20406

Chemical Properties

Cas. No. 214475-53-3

Formula $C_{22}H_{18}N_2O_4$ M.Wt 374.4

Solubility Soluble in DMSO Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

This linker which is stable to acid and bases is cleaved with high specificity under mild oxidative conditions in the presence of base and a nucleophile, yielding a range of C-terminal functionalities. The linker may be used to synthesize fully or part

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA