
Product Data Sheet

Product Name: 4-Methyl-benzhydrylamine resin (200-400 mesh, 0.70-1.30 mmol/g) .
HCl

Cat. No.: GA20422

Chemical Properties

Cas. No.

Formula

M.Wt

Solubility Soluble in DMSO

Storage

Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

Background

BHA resin turned out to be not sufficiently acid-labile for obtaining peptide amides with sterically hindered C-termini (e.g. Phe, Leu) applying standard cleavage conditions. Introduction of a 4-methyl substituent sufficiently increased the acid-labili

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA