
Product Data Sheet

Product Name: 5-FAM-HIV-1 tat Protein (47-57)

Cat. No.: GA20436

Chemical Properties

Cas. No. 1676104-81-6

Formula $C_{85}H_{128}N_{32}O_{20}$

M.Wt

1918.16

Solubility Soluble in DMSO

Storage

Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

FAM-HIV-1 Tat 47-57 was used as a fluorescent probe for studying the interaction between the p53 tetramerization domain and HIV-1 Tat protein.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA