
Product Data Sheet

Product Name: 5-TAMRA-Amyloid β -Protein (1-40)

Cat. No.: GA20440

Chemical Properties

Cas. No. 1802087-81-5

Formula $C_{219}H_{315}N_{55}O_{62}S$

M.Wt 4742.31

Solubility Soluble in DMSO

Storage Store at $-20^{\circ}C$

General tips For obtaining a higher solubility , please warm the tube at $37^{\circ}C$ and shake it in the ultrasonic bath for a while. Stock solution can be stored below $-20^{\circ}C$ for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

Anderson and Webb could verify using transmission electron microscopy that N-terminal labeling of A β 40 with TAMRA and other fluorescent dyes does not prevent the formation of protofibrils and amyloid fibrils of various widths.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA