
Product Data Sheet

Product Name: Ac-Phe-Lys-OH

Cat. No.: GA20634

Chemical Properties

Cas. No. 14287-21-9

Formula $C_{17}H_{25}N_3O_4$ M.Wt 335.4

Solubility Soluble in DMSO Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

Inverted aspartame-type sweetener. Adrover et al. used Ac-Phe-Lys as a model compound in their studies on protein glycation.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA