

---

**Product Data Sheet**

---

Product Name: Alarin (rat)

Cat. No.: GA20698

**Chemical Properties**

Cas. No. 909409-88-7

Formula  $C_{119}H_{199}N_{45}O_{35}$ 

M.Wt

2820.17

Solubility Soluble in DMSO

Storage

Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

Peptide belonging to the galanin family of peptides. Alarin exhibits potent dose-dependent vasoconstrictor and anti-edema activity in the cutaneous microvasculature. Alarin stimulated food intake and gonadotrophin release in male rats. glpbio provides t

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA