
Product Data Sheet

Product Name: Amyloid P Component (27-38) amide

Cat. No.: GA20716

Chemical Properties

Cas. No. 180387-75-1

Formula $C_{68}H_{107}N_{19}O_{17}S$

M.Wt

1494.78

Solubility Soluble in DMSO

Storage

Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

A wide variety of cells attach to the surface of polystyrene plastic dishes coated with this peptide, which corresponds to a fragment of the amyloid P component, a glycoprotein found in all types of amyloid deposits.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA