

---

**Product Data Sheet**

---

Product Name: Angiotensin II Antipeptide

Cat. No.: GA20764

**Chemical Properties**

Cas. No. 121379-63-3

Formula  $C_{42}H_{62}N_{10}O_{12}$  M.Wt 899.01

Solubility Soluble in DMSO Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

The sequence of this "anti-angiotensin II" peptide has been deduced from the antisense mRNA complementary to the human angiotensin II mRNA. The resulting octapeptide shows 50% sequence homology with angiotensin II. It acts, however, as an ant

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA