
Product Data Sheet

Product Name: APL1 β 28
Cat. No.: GA20778

Chemical Properties

Cas. No. 1233876-42-0

Formula C₁₀₉H₁₈₅N₃₁O₃₉S M.Wt 2585.92

Solubility Soluble in DMSO Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

Background

A highly sensitive biomarker for the production of A β 1-42 derived from the β APP-like protein APLP1. The peptide is generated by the same proteolytic mechanism as A β 42. Non-amyloidogenic APL1 β 28 can be detected in the CSF and its le

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA