
Product Data Sheet

Product Name: Atrial Natriuretic Factor (4-28) (human, bovine, porcine)

Cat. No.: GA20785

Chemical Properties

Cas. No. 95896-08-5

Formula $C_{112}H_{175}N_{39}O_{35}S_3$

M.Wt 2724.06

Solubility Soluble in DMSO

Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

Background

Antibodies raised to rat ANF (4-28) also reacted equally well with hANF (4-28). Anaritide was evaluated in the management of acute renal failure.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA