
Product Data Sheet

Product Name: Glpbio espresso mugs, pack of 6 assorted colors

Cat. No.: GA20792

Chemical Properties

Cas. No.

Formula

M.Wt

Solubility Soluble in DMSO

Storage

Store at -20°C

General tips

For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

Pack of 6 espresso mugs, 100ml/ 3.5 fl. oz. blue "B leading", cyan "B inspiring", green "B reliable", yellow "B brave", orange "B successful", magenta "B innovative" This merchandise is available

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA