
Product Data Sheet

Product Name: Biotinyl-(Gln¹)-LHRH

Cat. No.: GA20816

Chemical Properties

Cas. No. 218433-98-8

Formula C₆₅H₉₂N₂₀O₁₅S M.Wt 1425.64

Solubility Soluble in DMSO Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

This biotinylated LHRH derivative has been used in a functional analysis of the LHRH receptor ligand binding.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA