
Product Data Sheet

Product Name: Boc-Ala-Ala-Pro-Ala-pNA

Cat. No.: GA20883

Chemical Properties

Cas. No. 121570-42-1

Formula $C_{25}H_{36}N_6O_8$ M.Wt 548.6

Solubility Soluble in DMSO Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

Boc-AAPA-pNA, substrate for human neutrophil and hog pancreatic elastases.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA